

A preserved, isotonic, pH buffered, non-spermicidal, non-irritating, water-soluble lubricating jelly for all types of non-sterile medical procedures.

- Supports sperm cell functionality due to its isotonic (neutral) formulation
- Other products may be hypertonic, leading to loss of cellular fluid
- In a controlled study, fresh extended semen was mixed with CLARITY® A.I.
 Lubricating Jelly at 1% and 5% concentrations* with:
 - No effect on motility
 - No effect on viability (intact acrosome)

*Study final report available by request (multiple types of extenders were used).



AVAILABILITY

1 Gallon - Reorder No: 59008 12 oz. Bottle - Reorder No: 59007

www.aurorapharmaceutical.com

CLARITY is a Registered Trademark of Aurora Pharmaceutical, Inc.

FOR USE IN ANIMALS ONLY FOR EXTERNAL USE ONLY

For questions, please call 888-215-1256 Aurora Pharmaceutical, Inc. Northfield, MN 55057



University Research Highlights

A comprehensive Texas A&M University¹ study looking at five leading lubricating jellies (HR® Lubricating Jelly, K-Y™ Jelly, Therio-gel®, Priority Care® and CLARITY®) and the results highlight the importance of using caution when selecting an Al lubricant, even if lubricants are marketed as safe for this purpose.

- Exposure to 5% K-Y™, Priority Care®, or Therio-gel® lubricants yielded lower sperm quality
- The effect was most profound in K-Y™
- Most sperm quality measures were unaffected by different concentrations (1% or 5%) of HR® and CLARITY® lubricants with values similar to control
- Lubricant Therio-gel® tended to yield lower values for sperm lipid peroxidation; however,
 Therio-gel® increased sperm susceptibility to oxidative injury
- CLARITY® A.I. Lubricating Jelly yielded high values for all 5 endpoints, whereas H-R® yielded high values for only 4 endpoints; therefore, these lubricants might be the safest among those tested