

# CLARITY<sup>®</sup>



## Ultrasound Gel

Medium Viscosity Conductive, Salt & Alcohol Free, Water-Soluble Contact Medium

### INDICATIONS

CLARITY<sup>®</sup> Ultrasound Gel was developed as a coupling medium to facilitate transmission of the ultrasound energy from the machine head to the tissues.

- Highly conductive, clear (bubble free) transmission
- Long lasting, low absorption - no need to reapply as often
- Smooth consistency spreads easily without running or drying prematurely
- Non-greasy, non-staining, water-soluble clear formulation
- Contains no salt, alcohol or fragrance
- Non-corrosive - will not damage equipment

### APPLICATION

Ultrasound gel is used to avoid intense sound reflections caused by air pockets at the borderline between the ultrasonic probe and skin. The best ultrasonic media is free of air bubbles to ensure perfect sound transmission. Liquids like water and alcohols are suitable media, however, they are inappropriate with regard to their volatility and low viscosity. CLARITY<sup>®</sup> Ultrasound Gel was developed with fluidity and viscosity in mind.



### AVAILABILITY

1 Gallon - Reorder No: 59008  
16 oz. Bottle - Reorder No: 59007

### DIRECTIONS

Apply gel to treatment area. After use, remove with tissue or towel. Clean equipment after each use.

[www.aurorapharmaceutical.com](http://www.aurorapharmaceutical.com)  
CLARITY is a Registered Trademark of Aurora Pharmaceutical, Inc.

**FOR USE IN ANIMALS ONLY**  
**FOR EXTERNAL USE ONLY**



# CLARITY®

Ultrasound Gel

Medium Viscosity  
Conductive, Salt & Alcohol Free,  
Water-Soluble Contact Medium

For Animal Use Only

For External Use Only

1 Gallon (3.79 L)

## **DIRECTIONS:**

Apply gel to treatment area. After use, remove with tissue or towel. Clean equipment after each use.

Store at 15°– 30° C  
(59°– 86° F).

TAKE TIME  OBSERVE LABEL DIRECTIONS



MANUFACTURED IN THE USA  
ss000036 04/2020

For questions, please call 888-215-1256  
Aurora Pharmaceutical, Inc.  
Northfield, MN 55057